

## COMPARATIVE MATRIX FOR EXTRAORAL RANGE



Modalities and Features	CS 8100	CS 8100SC	CS 8100 3D	CS 8100SC 3D	CS 9300 Select	CS 9300C Select	CS 9300	CS 9300C
2D panoramic, TMJ, sinus	✓	✓	✓	✓	✓	✓	✓	✓
2D cephalometric	Upgradeable to ceph	✓	Upgradeable to ceph	✓	Upgradeable to ceph	✓	Upgradeable to ceph	✓
3D imaging	Upgradeable to 3D	Upgradeable to 3D	✓	✓	✓	✓	✓	✓
3D fields of view size (cm)	N/A	N/A	4x4, 5x5, 8x5, 8x9	4x4, 5x5, 8x5, 8x9	Selectable: 5x5, 8x8, 10x5, 10x10	Selectable: 5x5, 8x8, 10x5, 10x10	Selectable: 5x5, 8x8, 10x5, 10x10, 17x6, 17x11, 17x13.5	Selectable: 5x5, 8x8, 10x5, 10x10, 17x6, 17x11, 17x13.5

Applications	CS 8100	CS 8100SC	CS 8100 3D	CS 8100SC 3D	CS 9300 Select	CS 9300C Select	CS 9300	CS 9300C
General practice	•••	•••	•••	•••	••	••	•	•
Endodontics			•••	•••	••	••	•	•
Implantology	•	•	•••	•••	•••	•••	•••	•••
Periodontics	•	•	••	•••	•••	•••	•••	•••
Oral Surgery	•	•	•	••	•••	•••	•••	•••
Orthodontics		•••		••		•••		•••
Radiology center					••	••	•••	•••
ENT practice							•••	

• Good •• Better ••• Best

Please note: The bullet ratings listed above are simply recommendations based on feedback from our thought leaders and users. One bullet indicates the unit meets basic requirements for the application; three bullets indicate it is the best choice.