

Kodak dental films with a ClinAsept barrier

KODAK INSIGHT Dental Film with CLINASEPT Barrier

KODAK ULTRA-SPEED Dental Film with CLINASEPT Barrier

Kodak dental films with ClinAsept barriers are the only intraoral films with a built-in barrier to help manage infection control. Each intraoral packet is encased in a clear barrier that guards the film packet against contaminants that may be present in oral fluids. Use these products to minimize the risk of cross-contamination in your office.



Minimize cross-contamination and maximize confidence.

For real time savings and peace of mind, use Kodak intraoral dental films with ClinAsept barriers. Available in single- or double-film format, sizes 0 and 2 each come safely prepackaged with its own barrier. You don't have to do anything extra before taking the radiograph, except feel confident that the risk of cross-contamination in the darkroom is minimized.

Get in touch today

Make Kodak intraoral dental films with ClinAsept barriers part of your key infection control practices. For speed and convenience, they're a real plus in your dental practice.

How to use prepackaged films with barriers



Ordering Information

For KODAK Dental Films with a CLINASEPT Barrier

| Code | Size | Quantity Per Package | CAT Number |
|--|--------|----------------------|------------|
| Kodak InSight dental films with a ClinAsept barrier | | | |
| IP-21C | Size 2 | 100 1-film packets | 153 9931 |
| IP-22C | Size 2 | 100 2-film packets | 156 0390 |
| IP-01C | Size 0 | 75 1-film packets | 120 0328 |
| IP-02C | Size 0 | 75 1-film packets | 124 7410 |
| Kodak Ultra-speed dental film with a ClinAsept barrier | | | |
| DF-58C | Size 2 | 100 1-film packets | 171 7131 |
| DF-57C | Size 2 | 100 2-film packets | 898 7877 |
| DF-54C | Size 0 | 75 1-film packets | 144 5360 |

For more information, call **1-800-933-8031** or visit: www.kodakdental.com



Carestream Health, Inc.

© Carestream Health, Inc., 2008.
InSight, Ultra-speed, Super Poly-Soft and ClinAsept are trademarks of Carestream Health, Inc.
The Kodak trademark and trade dress are used under license from Kodak.

D3-37 883 7668
2009-07-22