

Carestream Dental strongly adheres to compliance with Microsoft Windows operating systems and hardware compatibility. Only operating systems and hardware listed in this document have been fully tested and certified on WinOMS CS. Hardware and software listed in this table is recommended for new purchases. Before purchasing any equipment, please review the latest system requirements located online at: <http://www.carestreamdental.com/documentation/sysReqsWinOMS.pdf>.

<p align="center">Dedicated Server</p> <p align="center">(reserved to run the DATABASE SOFTWARE)</p>	<p>Hardware</p> <ul style="list-style-type: none"> ▪ X86 or x64 CPU Dual-Core ▪ 4 GB RAM (8 GB Recommended) ▪ 100 GB available disk space ▪ 3D accelerated video card ▪ DVD drive ▪ Ethernet 10/100/1000 network card ▪ 2 or more powered USB 2.0 ports ▪ Additional PCI Express, AGP or USB 2.0 expansion slots may be required 	<p>Supported Operating Systems</p> <p>Windows Server 2008, SP 2, 32-bit and 64-bit Windows Server 2008 R2 Windows Server 2003 Standard Edition, SP 2, 32-bit</p>
<p align="center">Workstation</p>	<p>Hardware</p> <ul style="list-style-type: none"> ▪ X86 or x64 CPU Dual-Core ▪ 2 GB RAM (4 GB Recommended) ▪ 40 GB available disk space ▪ 3D accelerated video card ▪ DVD drive ▪ Ethernet 10/100/1000 network card ▪ 2 or more powered USB 2.0 ports ▪ Port(s) appropriate for printer and scanner connectivity ▪ Standard LCD monitor capable of 1280x1024 screen resolution. ▪ Additional PCI Express, AGP or USB 2.0 expansion slots may be required ▪ DB9 serial port for Anesthesia Vitals Monitor 	<p>Supported Operating Systems</p> <p>Windows 7 Professional, 32-bit and 64-bit Windows Vista Business, SP 2, 32-bit and 64-bit Windows XP Professional, SP 3, 32-bit</p>
<p align="center">Tablet PC</p> <p align="center">(intended for use with the Point-of-Care Module)</p>	<p>Hardware</p> <ul style="list-style-type: none"> ▪ Intel Pentium M 2.0 GHz or faster ▪ 1 GB RAM ▪ 40 GB available disk space ▪ DVD drive ▪ Ethernet 10/100/1000 network card ▪ Wireless WLAN card supporting 802.11b or faster ▪ Port(s) appropriate for printer and scanner connectivity ▪ Point-of-Care requires a minimum 768x1024 screen resolution 	<p>Supported Operating Systems</p> <p>XP Tablet PC edition, 2005, SP3 Windows Vista Business, SP 2, 64-bit Windows 7 Professional, 32-bit and 64-bit</p>
Additional and Add-on Modules for WINOMS CS Software		
<p align="center">Anesthesia Records Module</p>	<p>CUSP Compatible Criticare 8100E nGenuity . (Contact Criticare customer service to determine whether your vitals monitor is CUSP compatible), 9-pin Male/Female Null-modem Serial cable. Criticare DB-9 Bluetooth device is supported. A USB adaptor is not supported.</p>	
<p align="center">Data Export</p>	<p>Microsoft Excel 2007, SP2 or 2010 alone or as a part of a Microsoft Office 2007 or 2010 Edition. Having multiple versions of Microsoft office throughout the practice may cause document incompatibilities and is not recommended .</p>	
<p align="center">Digital Forms Module</p>	<p>Microsoft Word 2007, SP2 or 2010 (32bit ONLY) alone or as a part of a Microsoft Office 2007 or 2010 (32bit ONLY) Edition. Having multiple versions of Microsoft office throughout the practice may cause document incompatibilities and is not recommended. Topaz Systems SignatureGem LCD 4x3 model: T-LBK755-HSB or T_LBK755-HSB-R. Topaz Systems SignatureGem LCD 1x5 Model: T-LBK462-HSB.</p>	
<p align="center">Image, Document, and ID Card Scanning</p>	<p>Scanners listed on page 2 have been certified with WinOMS CS software. Any TWAIN compliant scanner is acceptable. Network scanners are not supported. Duplex and multi-page features supported for documents and ID cards. Image scanning requires a flatbed scanner with a full-sized transparency adapter.</p>	
<p align="center">Optical Ticket Scanning</p>	<p>Epson GT 1500 Workforce, or Epson Perfection v700 Photo</p>	
<p align="center">Printers</p>	<p>Local and network printing are supported. Printers listed on page 2 have been certified with WinOMS CS software, but any HP LaserJet printer using the HP PCL 6 driver is acceptable. For chart labels choose a Laser or Dot Matrix Label printer.</p>	
<p align="center">eClaims, Automated Statements</p>	<p>High-speed broadband internet connection or US Robotics 56K EXT modem</p>	
<p align="center">ePayments</p>	<p>High-speed broadband internet connection (recommended) and Internet Explorer 7 or US Robotics 56K EXT modem. MAGTEK credit card reader (USB or serial) only available from Carestream Dental.</p>	
<p align="center">eVerifications</p>	<p>Insurance Eligibility and Claim Status Services are supported. High-speed broadband internet connection and Internet Explorer 6 or higher are required.</p>	
<p align="center">KODAK OMS Imaging Software</p>	<p>Image Server, Image Capture Station Refer to this location for detailed specifications: http://www.carestreamdental.com/documentation/sysReqs/OMS_Imaging.pdf</p>	
<p align="center">Mail Merge</p>	<p>Microsoft Word 2007, SP2 or 2010 or as a part of a Microsoft Office 2007 or 2010 Edition. . Having multiple versions of Microsoft office throughout the practice may cause document incompatibilities and is not recommended.</p>	
<p align="center">MasterWorks LearnCenter</p>	<p>High-speed broadband internet connection, Internet Explorer 7.0 or higher. MasterWorks LearnCenter additionally requires a sound card, speakers, headset with microphone, Adobe Flash Player 9, Java/JRE 5.0</p>	
<p align="center">NEA FastAttach</p>	<p>NEA FastAttach software version 3.2 or higher. KODAK E-Claim Service. Windows XP Professional or Windows Vista without Terminal Services or Citrix. High-speed broadband internet connection, Internet Explorer 6.0 or higher.</p>	
<p align="center">Point-of-Care Module</p>	<p>Tablet PC, Wireless LAN</p>	
<p align="center">PEARL</p>	<p>Supported devices are the iPhone 3G and 3Gs, and BlackBerry Bold, Torch and Curve. iPod Touch (current generation) devices will run PEARL, but require a Wi-Fi data connection and have no native phone/text support.</p>	
<p align="center">QuickBooks Software Module</p>	<p>QuickBooks Pro 2010.</p>	
<p align="center">TRUFORM Software Module</p>	<p>Adobe Acrobat version 9.0 or higher (full version), Topaz signature capture plug-ins for Adobe Acrobat, Internet connection. Topaz Systems SignatureGem LCD 4x3 model T-LBK755-HSB or T_LBK755-HSB-R. SignatureGem LCD 1x5 Model T-LBK462-HSB</p>	

The following additional hardware and software has been certified with previous versions of WinOMS CS software:

Operating Systems:

Microsoft Windows XP Professional, SP2, 32-bit
Microsoft Windows Vista Home Premium, Enterprise, Ultimate 32-bit
Windows XP Tablet PC Edition, 2005
Microsoft Windows 7 Professional, Ultimate 32-bit

General Printers:

HP Color LaserJet CP4005 Printer Series with HP PCL 6 driver
HP LaserJet P4515 or 3005d Printer Series with HP PCL 6 driver
These printers, while discontinued by HP, are approved for use with WinOMS CS: 1012,1300,1320,2200,2300,2430,4000,4050,4100, 4200,4250,4350.

Label Printers:

Brother QL-500
Seiko SLP 420, 430, 440, 450

TRUFORM Software Module:

Adobe Acrobat 6 , 7 , 8.1 (full version)

Image Scanners (with Transparency Adapters):

Epson Perfection v700 Photo
These scanners, while discontinued by their manufacturers, are approved for use with WinOMS CS: HP ScanJet 4 C/T; Microtek Scanmaker II ; Epson Expression 800, Epson Expression 1600, Epson Expression 1680 Professional, Epson Perfection 4990 Photo

Anesthesia Vitals Monitors:

Criticare Connect Bluetooth Adapter

Criticare 8100H (nCompass)
Criticare 8100E (nGenuity)
Criticare 506 NT3 (not ECG compatible)
These monitors, while **discontinued** by Criticare, are approved for use with WinOMS CS: Criticare Poet Plus 8100, Criticare 507EL Scholar (not ECG compatible)

Digital Cameras:

KODAK P712 Dental Digital Photography System
KODAK 6490 or 7590 Systems
Canon EOS Rebel XSi

Document and ID Card Scanners:

Epson WorkForce GT-1500 Scanner
Fujitsu fi-6130
HP ScanJet7800
Kodak SCANMATE i1120
Ambir DS487-AS
CSSN Scanshell 800DXN
Epson Perfection v700 Photo
Visioneer Strobe XP 300
These scanners, while discontinued by Epson, are approved for use with WinOMS CS: Expression 1600, Expression 1680 Professional, Perfection 4990 Photo.
These scanners, while **discontinued** by HP, are approved for use with WinOMS CS: 5200c,6200c,6300c,7400c.
Visioneer Strobe XP 200, while discontinued, is approved for use with WinOMS CS.

Mail Merge Software:

Microsoft Word XP, 2003

Data Export Software:

Microsoft Excel XP, 2003

Accounting Software:

QuickBooks Pro 2008, 2009, 2010

Internet Browser:

Internet Explorer 6, 7, 8,9

Anti-Virus Software

Internet security software providing anti-virus and other malware protection for your WinOMS CS server and workstations is recommended. Compare products for system performance.

Notes & Recommendations:

- **Microsoft Office:** In order to use the **digital forms module**, users must have one of the following: Microsoft Word 2007, SP2 or 2010 (32-bit version only) or as a part of a Microsoft Office 2007 or 2010 (**32bit ONLY**) Edition (Standard, Small Business or Professional editions, 32-bit version only), Microsoft Office 2003 Professional or Microsoft Word 2003 Professional. Other Microsoft Office versions approved for **Mail Merge ONLY** with WinOMS CS software are: Microsoft Office 2003 (Standard, Basic and Small Business), Office 2002 (XP), 2000, or 97. **WinOMS CS does not support Home, Student editions of Microsoft Word or the Office Suite of products, or Office 2010 64-bit versions.** Having multiple versions of Microsoft office running in the same office location is not recommended due to document incompatibilities.
- **Insurance Forms:** Contact RR Donnelley for insurance related forms: 800 401 9771
- **QuickBooks:** Standard, Premier, Professional and Enterprise software editions are the only versions supported by WinOMS CS software. The Simple Start, Online and Canadian software versions are not approved for use with WinOMS CS software.
- **Backups and Uninterruptible Power Supplies:** Carestream Dental LLC suggests that you perform a nightly backup on removable media and keep that media in a fireproof safe off-site. The database engine must be stopped to make a proper copy. Carestream Dental LLC also recommends that you have a UPS with shutdown software installed with the server machine.
- **Operating System Updates:** Make sure that all servers and workstations are Microsoft Certified for their respective operating systems. Each operating system should have the latest service packs and critical updates installed.
- **Hardware and Software Performance:** System Requirements are updated regularly. Hardware upgrades (processor, disk, memory, graphics and/or network) may be required to take full advantage of the new features in the latest versions of the software. In general, performance can be improved by doing one, some or all of the following: Increasing the amount of memory (RAM) on a workstation or server; replacing the hard drive with a newer, larger, faster model; upgrading the CPU to a faster speed. Results are not guaranteed.
- **Printer Performance:** To improve printer response on Windows XP Professional, turn off Windows Advanced Text Services at Control Panel>Regional and Language Options>Languages tab>Details>Advanced tab.
- **Server:** Carestream Dental LLC strongly recommends running the WinOMS CS database on a dedicated server.
- **Internet Connection:** Although it is not necessary for WinOMS CS to have an internet connection to fully-function, Carestream Dental LLC strongly recommends that there be an always-on internet connection available to your network. This allows the Carestream Dental Technical Solutions team to provide assistance quickly and easily to your software and hardware when needed. If this kind of internet connection is not feasible for your practice, having a modem and remote control software (such as PCAnywhere) installed on your server will suffice.
- **Microsoft Windows Vista:** WinOMS CS versions 7.0.18 and greater support the Microsoft Windows Vista Operating System. Please check with your third-party hardware and software vendor to ensure that their products are compatible with Windows Vista before installing that operating system. Microsoft Word or Microsoft Office 2007 Standard, Professional, Small Business or Enterprise Editions are the only compatible versions for use with Windows Vista operating systems.
- **Microsoft Windows 7:** WinOMS CS versions 8.1.0 and higher are compatible with Windows 7, all editions except Home Premium. WinOMS CS is only compatible with Microsoft Office 2007 running on Windows 7. Please check the Kodak OMS Imaging Software System Requirements for Windows 7 compatibility before choosing to run WinOMS CS integration with Kodak OMS Imaging software on Windows 7. Please check with your E-services or third-party hardware and software vendor to ensure that products are compatible with Windows 7.
- **Anesthesia Record Integration:** All Criticare Vital Sign Monitors contain a Female serial port connection while all standard workstations contain a Male Serial Port connection. Purchasing a 9-Pin **Male/Female** NULL-Modem serial cable is required. A USB adaptor is not supported.
- **Document and ID Card Scanning:** Automatic document feed (ADF) is recommended. Download and install the latest TWAIN driver from the manufacturer's website before using. Apply the latest driver for Ambir DS487-AS from www.ambirtechnology.com/drivers/487patch.exe.
- **General Printing:** Download and install the latest HP PCL 6 driver from the manufacturer's website before using.
- **Patient Workspace:** **Patient Workspace requires a minimum 1024 x 768 screen resolution. For best visibility of Patient Workspace and Patient Dock**, select monitors and video cards that deliver 1280 x 1024 resolution. To run Patient Workspace on a tablet PC in portrait orientation, a minimum 1024 horizontal screen resolution is required.
- **Topaz Signature Pads:** "BHSB-R"= **Both** Serial and USB "BHSB-R"= RoHS compliant or lead free