

Kodak

Dental Imaging Software

Designed to work with our family of imaging devices, Kodak dental imaging software helps you manage and analyze all of your digital dental images. Each patient's images are well-organized for easy access and can be quickly analyzed to help you provide your patients with the highest level of professional care.



Easy to use and integrate.

Simple, efficient, and intuitive, Kodak dental imaging software provides tools to help you make fast, confident diagnoses. It can be used with its own patient file as a stand alone program, or it may be integrated with your existing practice management software.

The right tools, helping you work smart to enhance patient care.

Work smarter and more efficiently with image management and processing tools designed specifically for the Kodak RVG 6100 digital radiography system, the Kodak 8000 digital panoramic system and the Kodak 1500 intraoral camera.

- The image processing functions feature digital x-ray tools and controls such as contrast and brightness, three pre-defined optimized image viewing modes (periodontics, endodontics, and DEJ), high-resolution filter, highlight, zoom, measurements, negative video, pseudo-colors, sepia filter, pseudo 3D, relief filter, density analysis tools, drawings and comments, and the optional Kodak cosmetic imaging software.
- In the display mode, choose between viewing your individual images full screen or actual image size. View and compare your digital x-rays in operative radiology mode, four image bitewing mode, or within a variety of full mouth series mounts.
- Keep a detailed image log per patient, conduct a detailed image search and delete, archive, and import/export images.
- With the printing and transfer function, print all image types, reports, and send images by email (with or without our image viewer).

Improved communication and diagnoses.

Clean, clear images and powerful image manipulation tools can aid in more consistent patient communication and assist in more efficient and concise patient diagnoses. Combined with the Logicon® caries detector, to help detect the presence and depth of penetration of proximal caries on permanent teeth, Kodak dental imaging software can help enhance all of your digital imaging capabilities.

Support and guidance you can count on.

We're committed to providing the ongoing guidance and support you need to ensure a smooth transition to our products.

Call 800-944-6365 today to take advantage of the practice enhancing tools.

FOR MORE INFORMATION, VISIT:
www.carestreamdental.com



Carestream Dental

©Carestream Health, Inc., 2010.
The Kodak trademark and trade dress are used under license from Kodak.
RVG is a trademark of Carestream Health.

2530 10/15/10

