

# WHAT'S NEW IN CS IMAGING VERSION 8



## 1 POWERFUL DASHBOARD IMAGE GALLERY

- Gives an overview of the patient history, including thumbnails of 2D images, 3D images and CAD/CAM data
- View multiple display and filter modes
- Preview large-size images, including 3D images and CAD/CAM data that can be rotated
- Attach documents, images and videos to patient records (Word files, PDF, JPEG, TIF, MPEG)
- Access acquisition interface, add comments, share, print and export images directly from the Dashboard

### Benefits

Find images without opening the review workspace or any 3D software drastically improves practice workflow.



## 2 FASTER IMAGE REVIEW

- New auto-arrange function automatically organizes and resizes images for optimal display without image overlap

### Benefits

Faster, easier image review with fewer clicks allows users to save time, and is ideal for multiple images review (FMS, bitewing, endodontic treatments and more).



## 3 INTEGRATION OF 3D IMAGING AND CAD/CAM

- Display any 2D images, 3D images and CAD/CAM data simultaneously without switching between programs
- Automatically store extracted CBCT slices, CBCT screenshots and other screenshots generated in our 3D and CAD/CAM software in the patient history
- Display multiple CBCT volumes side-by-side for basic image comparison

### Benefits

Manage all data within one software. Simplifies treatment planning, case presentation and pre-/post-op monitoring streamlines the clinical workflow.



## 4 ENHANCED INTEGRATION AND INSTALLATION

- Software can be installed without drivers
- No change to third-party practice management software integration
- Low system requirements (similar to CS Imaging version 7)

### Benefits

Minimizes installation time thanks to easy integration. No need to upgrade existing computer.



## 5 MODERN, USER-FRIENDLY DESIGN

- Tools are easier to access and images always well-organized on the screen
- Black user interface enhances the diagnostic experience
- Extended use of right-click commands
- Invisible image frame windows leave more space for reviewing images

### Benefits

Improved user experience and chairside use streamlines the workflow and saves time.



## 6 ENHANCE MULTIPLE IMAGES AT ONCE

- Enhance multiple images simultaneously using the same processing on all selected images (FMX, bitewing, endodontic exams and more)

### Benefits

Save time when diagnosing a series of X-rays.



## 7 IMPROVED IMAGE SHARING

- Import and export image via simple drag and drop
- Export images directly into an email, document or folder

### Benefits

Facilitates case collaboration and data sharing.



## 8 NEW DARKROOM MODE

- Displays images on an uncluttered full screen
- Maximizes the image area while keeping all tools and image acquisition accessible

### Benefits

Improved image review and case presentation aids in diagnosis, treatment planning and patient communication.



## 9 IMPROVED MEASUREMENTS, DRAWINGS AND ANNOTATIONS TOOLS

- Make measurements, drawings and annotations directly on images
- Compare measurements made on multiple selected images
- Redo / undo changes on measurements and drawings
- Add comments directly on images and 3D data in the Dashboard
- Trace cephalometric images in less than 10 seconds, five times faster than before

### Benefits

Faster, more intuitive measurements, drawings and comments for better recordkeeping and treatment planning.



## 10 EFFORTLESS RECORD CORRECTION

- Re-assign images acquired incorrectly under the wrong patient name
- Move images to another patient in a single click

### Benefits

Ensure accuracy of records by easily moving images acquired in the wrong patient file to the correct chart.