



## CS 2200 INTRAORAL X-RAY SYSTEM

# Performance. Operational Ease. Safety.

**Obtain high-quality images consistently,  
while improving safety**

Compact, intuitive and easy to use, the CS 2200 offers the superior performance you need while providing improved safety for your patients and staff—year after year. The unit's high frequency DC X-ray generator ensures delivery of the image quality you're looking for, while the compact, lightweight design accommodates the space requirements of almost any practice.



800.944.6365  
[carestreamdental.com](http://carestreamdental.com)

# CS 2200: Exceeding Your Expectations for 2D Imaging

## MAXIMUM USABILITY

The intuitive, handheld timer allows users to quickly select correct settings or choose from pre-programmed options

## MINIMAL RADIATION EXPOSURE

The CS 2200 system automatically determines the optimal dose, maximizes the amount of useful X-rays and generates up to 30% less radiation\*

## IMAGE OPTIMIZATION

Equipped with adjustable tension (60 or 70 kV), the CS 2200 system lets you accentuate contrast and levels of gray shading, according to your diagnostic needs

## HIGH IMAGE QUALITY

With a focal spot of 0.7 mm, the system produces sharp, detailed images for diagnostic excellence



## IMAGING COMPATIBILITY

Compatible with both film and digital imaging applications, the CS 2200 system works the way you do

## CONFIGURATION VERSATILITY

With multiple installation configurations available, the CS 2200 system easily adapts to your practice's needs\*\*

## PRECISE POSITIONING

The stable arm and lightweight tube head ensure correct, effortless alignment and reduce the risk of cone-cutting images

We are the only dental company to engineer and manufacture X-ray generators from start to finish, enabling us to maintain complete control over the standards and quality of each system.

\* When compared to conventional X-ray generators

\*\*The CS 2200 system's wall mount was specifically designed to facilitate the replacement of the former Trophy Irix unit to minimize the impact and cost to your practice.