



CS 7200 IMAGING PLATE SYSTEM

High Image Quality. Workflow Simplicity.

Step into digital with the ideal chairside companion.

Designed for maximum simplicity, the CS 7200 imaging plate system lets you expose, scan and view. Images open on your computer screen, and the plate automatically erases when you're finished. You can capture images, scan plates and review images—all from the dental chair.

The CS 7200 produces impressive image quality (19 lp/mm) and includes advanced image processing tools that optimize contrast. CS Adapt's preprogrammed filters enable you to enhance images according to your diagnostic needs. Advanced algorithms generate images with details that are more visible and helpful in the diagnostic process.



800.944.6365
carestreamdental.com

Designed to fit your workspace— and budget

With its affordable price, high reliability and low maintenance costs, the CS 7200 is a cost-effective solution for any practice. Because it's digital, you can immediately eliminate your film and chemical expenses. The plates can be reused hundreds of times, and replacements—when needed—are priced reasonably. The space-saving design and quiet scanning process make the system ideal for chairside use.

GREATER PATIENT COMFORT

Increase comfort for your patients with a quiet scanning process and thin, flexible plates that are easy to position.

COMPLEXITY-FREE IMAGING

Acquire images directly in the imaging software; reuse plates hundreds of times.

MORE PERSUASIVE CASE PRESENTATION

Use CS Adapt and the imaging software to enhance your images and improve patients' understanding of treatment plans.

COVERAGE FOR DAILY INTRAORAL INDICATIONS

Accommodate your patients with plates for pediatric, periapical and bitewing exams in three sizes: 0, 1 and 2.

SIMPLE CONNECTIVITY

Connect effortlessly via USB.

SPACE-SAVING DESIGN

Enjoy a system that maximizes your workspace.

